



Mathematisch-Naturwissenschaftliche Fakultät Kepler Center for Astro and Particle Physics



Cavity enhanced long light path attenuation length measurement

An alternative approach in measuring the attenuation in the liquid scintillator for JUNO

Tobias Heinz, 18.09.17





Forschergruppe JUNO



Outline

- Motivation
- Properties of optical cavities
- Measurement of attenuation length
- First results
- Next steps



Motivation

- Attenuation length cannot be measured precisely with commercial available UV-vis spectrometers due to small decrease in intensity (the intensity drops to 99,6% of initial intensity after 10 cm in medium with attenuation length of 25 m)
- Need long light path to be able to measure attenuation length of LS
- But: very impractical to build cells with lengths of several meters
- Solution: use optical cavity to extend light path



Optical Cavity

An optical cavity is an arrangement of mirrors so that light is reflected multiple times

- An optical cavity is characterized by mirror reflectivities R_1, R_2 , the radii of curvature r_1, r_2 and the distance between both mirrors L
- The light path is only closed inside of cavity for certain combinations of r_1, r_2 and L

Stability criterion

$$0 \leq \left(1 - \frac{L}{r_1}\right) \left(1 - \frac{L}{r_2}\right) \leq 1$$



Optical Cavity

- Photons are reflected multiple times inside of cavity
 - Depending on the mirror reflectivities R_1 , R_2 the average light path is extended by a certain factor







LPE factor for cavity with effective mirror reflectivity $\sqrt{R_1R_2}$



Optical Cavity

- Photons are reflected multiple times inside of cavity
 - Depending on the mirror reflectivities R_1 , R_2 the average light path is extended by a certain factor

Light path extension factor
$$\epsilon = \frac{2R_1R_2}{1-R_1R_2}$$

• Effective average light path is given by:

Effective light path

$$L_{\rm eff} = (1 + \epsilon)L$$



Parameters of cavity mirrors

In Tübingen we use two identical concave cavity mirrors

- Radius of curvature: -50 cm
 - Possible cavity lengths from 0 to 1 meter
- Reflectivity: 0.95
 - LPE factor: ~18.5

Effective light path up to ~19.5 m possible



Measurement of attenuation length

The output intensity is a superposition of waves with various numbers of round trips

- Attenuation length cannot be measured with Beer's law

Two possibilities:

- Modulated intensity method:
 - initial intensity is modulated sinusoidally
 - output intensity is phase-shifted sine wave
 - Phase-shift occurs due to attenuation inside of cavity
- Ring down measurement:
 - build up intensity inside of cavity
 - turn off light source
 - output intensity will drop exponentially
 - time constant of exponential decrease depends on attenuation



CELLPAL measurement

Modulated intensity method was chosen:

- > ability to sweep modulation frequency
 - >more data due to relation between phase-shift and modulation frequency
- > the attenuation has an impact on amplitude and phase-shift of output intensity

>two ways to determine attenuation length with one measurement



Schematic illustration of setup





CELLPAL measurement

Input intensity

 $I_{\rm in}(t) = I_0 \left[1 + \alpha \cdot \sin\left(-\omega t\right)\right]$

measured with detector in front of cavity

Output intensity $I_{out}(t) = I'_0 [1 + \alpha' \sin(-\omega t + \phi)]$

measured with detector behind cavity



CELLPAL phase-shift measurement

Phase-shift between detector signals

$$\Delta \phi = rac{2\pi f \Delta L}{v_{g}} + \arctan\left(rac{A_{
m eff} \sin\left(rac{4\pi f L}{v_{g}}
ight)}{1 - A_{
m eff} \cos\left(rac{4\pi f L}{v_{g}}
ight)}$$

with
$$A_{\rm eff} = R_1 R_2 \exp\left(-\frac{2L}{L_{\rm att}}\right)$$

- *f*: modulation frequency
- ΔL : net distance between both detectors
- v_g : group velocity
- Latt: attenuation length



CELLPAL amplitude measurement

Output amplitude

$$l'_{0} = \frac{(1 - R_{1})(1 - R_{2})\exp\left(-\frac{L}{L_{att}}\right)}{1 - A_{eff}} l_{0}$$

$$\alpha' = \frac{(1 - A_{eff})\alpha}{\sqrt{1 - 2A_{eff}\cos\left(\frac{2\omega L}{V_{g}}\right) + A_{eff}^{2}}}$$
with $A_{eff} = R_{1}R_{2}\exp\left(-\frac{2L}{L_{att}}\right)$



• Light source: High Power LED (420 nm) coupled to fiber







setup placed inside darkbox with nickel-copper shielding











First results – Laser as light source/photodiode as detector





First results – Laser as light source/photodiode as detector

- The measured data points were fitted with the function for $\Delta \phi(f)$:
 - values used in calculation:
 - cavity length: 0.25 m
 - attenuation length in air: 5000 m
 - speed of light in air: 299 700 000 m/s
 - fit parameters:
 - mirror reflectivity
 - net distance between detectors

 $\begin{aligned} & \text{Results of fit} \\ & R = 0.956 \pm 0.004 \\ & \Delta L = 0.78\,\text{m} \pm 0.07\,\text{m} \end{aligned}$

Value of manufacturer $R = 0.95 \pm 0.005$

method seems to work!



First results – Amplitude measurement in Laser/PD setup





First results – Amplitude measurement in Laser/PD setup

- The measured data points were fitted with the function for $\frac{I'_0 \alpha'}{I_0 \alpha}(f)$ with the same values used in calculation as before
- Mirror reflectivity as fit parameter

Result of amplitude fit
$$R = 0.956 \pm 0.002$$



First results with current setup





First results with current setup

Same values used in calculation as before

LED/PMT setup result

 $R = 0.953 \pm 0.006$



First results with current setup

Same values used in calculation as before

LED/PMT setup result $R = 0.953 \pm 0.006$

Laser/PD setup result

 $R = 0.956 \pm 0.004$

Result of amplitude fit $R = 0.956 \pm 0.002$

Value of manufacturer $R = 0.95 \pm 0.005$

reproducable results with both setups!



Next steps

- Build cavity which is able to be filled with LS
- Reflectivity of mirrors in LS should be known precisely
 nood experimental method to measure mirror reflectivity
 - need experimental method to measure mirror reflectivity in liquid with refractive index >1



Thank you.

Contact:

Tobias Heinz

Kepler Center for Astro and Particle Physics University of Tübingen Auf der Morgenstelle 14 72076 Tübingen · Germany Phone: +49 7071 29-76280 tobias.heinz@uni-tuebingen.de







CELLPAL ring down time measurement

- After light source is turned off, the output intensity drops exponentially: $I(t) = I_0 e^{-t/\tau}$
- Ring down time depends on mirror reflectivity R, cavity length L and attenuation length L_{att}

Ring down time $\tau = \frac{L}{c\left[1 - R + \frac{L}{L_{att}}\right]}$



Frequency dependence of phase-shift for different values of attenuation length in 40cm cavity





Frequency dependence of amplitude for different values of attenuation length in 40cm cavity

