

ML Feedback for HI Jena Laser plasma accelerators


Alexander Sävert

Helmholtz Institute Jena

Jena, 14.12.21

High Energy Lasers

JETi200



Wavelength: Medium: Energy on target: Pulse duration: Peak power: Repetition rate: Probe beam: 800 nm Ti:Sapphire up to 5 J (4 J typical) 17 fs (20 fs typical) 300 TW 5 Hz 5 fs

POLARIS



Wavelength:1030 nmMedium: Yb:CaF2, Yb:GlassEnergy on target:up to 20 J (Comp. limited)Energy uncompressed:54 JPulse duration:>90 fsPeak power:>200 TWRepetition rate:1/40 Hz

JETi 200 – 200 TW Jenaer Titan Saphir Laser



Double CPA system with XPW filter installation in 2013-2015




Optimizing Laser driven particle acceleration



Laser parameter for optimization

- energy
- focus (shape, position)
- pulse duration (chirp)
- temporal intensity contrast



Experimental parameters for optimization

- target gas
- target length
- focus position
- pre plasma conditions


PIC simulations



www.hi-jena.de


4

Acceleration electron bunches to GeV level



Electron bunch spectrum


Cryo cooled power amplifier



- 100 W pumping power
 @532 nm
- target temperature 90 K to reduce the thermal lens of the Ti:Sapphire crystal
- Helium expander cryo head



Beam stabilization (short term drift)



- Short term jitter > 30 Hz
- Long term drift only depends on environment (temperature, humidity, air pressure)

HI JENA Helmholtz Institute Jena

Active beam stabilization



fast mirror ~ 1 kHz

For online measurement: Use transmitted light through high reflective mirrors and focus the beam on the 4Q diode.


Challenge:

- cw- pilot beam for continuous signal

- single shot full power beam for experiments
- → fast shutter to protect 4Q diode, but robust enough to withstand full power focused beam (closing time 30 ms, mirror stays put)

HI JENA Helmholtz Institute lena

Active beam stabilization with cryo cooling



Active beam stabilization

- needs an continuous online signal
- without pilot beam \rightarrow put in hold position
- Shutter closing time > typical jitter period (>30 Hz)
- → no active stabilization during high power shots




Methods for prediction

• has to be fast enough (analysis + mirror motion+ spare time) < 30 ms

Machine learning / Neural network	Autocorrelation search
"real" prediction	searching for similar signals from a database
	fast but not accurate
input layer hidden layer 1 hidden layer 2 output layer	




Results: neural network (vertical dev. medium jitter)



HI JENA Helmholtz Institute Jena


Results: neural network (horizontal dev. small jitter)





Eliminating the root cause

Remove the source of the vibrations ;)




+ use NN knowledge to reduce beam pointing jitter even further



Outlook for 2022

- So far: logging of laser parameter
- In 2022: Using Tango environment for control of experiments but a lot of legacy LABVIEW programs





HI Jena extension (finished in 2022)




September 2021

Target Area Fraunhofer (TAF)



Target Area Fraunhofer -TAF





High intensity lasers @ HI Jena (2022)

JETI ONE

JETi200

POLARIS







Wavelength: Energy on target: Pulse duration: Peak power: Repetition rate:

0.8-7 μm max. 7 mJ few cycle 1 kHz Wavelength: Energy on target: Pulse duration: Peak power: Repetition rate:

800 nm 5 J 17 fs 300 TW 5 Hz Wavelength: Energy on target: Pulse duration: Peak power: Repetition rate:

1030 nm 16 J (54 J) 100 fs 160 TW 1/50 Hz

Goal: synchronisation of all beams <20 fs (rms) and < 2 μrad (rms) in new target area



Target Area Fraunhofer -TAF



Targetarea Fraunhofer

HI JENA Helmholtz Institute Jena

Synchronisation & control: toddler steps (2022)



SLAVE

RF-synchronisation ca. 200 fs (rms)

active beam stabilization later

SLAVE HIJENA Helmholtz Institute Jena

Laser control system install base

- POLARIS control system adapted from PHELIX /GSI but most LABVIEW legacy
- JETi200 custom made CCM from Amplitude technologies
- JETIONE will be adapted

One control system needed!

HI JENA Helmholtz Institute Jena

Outlook beyond 2022



- Tango control system for controlling the laser, experiments and safety system (Tango + Siemens SPS)
- first steps: experiments & safety system (needed for new Target area,)
- Challenges -> Labview environment (Data aquisition, analysis, acceptance) keep GUI

