
Status of ECAL-P mechanics for Test Beam 2025

Grzegorz Grzelak, Piotr Zbińkowski, Filip Żarnecki

Faculty of Physics University of Warsaw




LUXE ECAL-P weekly zoom meeting 10 April 2025



• load > 50 kg, resolution \sim 0.5 [deg] (manual wrt the scale), better with laser pointer wrt the mirror \rightarrow work in progress

Test Beam configuration as discussed at DESY in Feb'25




 Sensors at positions (2,3) to be shifted to (3,4) at the center on rotation axis (accepted by Jakub)

Old FCAL Tungsten Plates ($14 \times 14 \text{ cm}^2$)




- 6 plates removed from permaglass frames and cleaned (Z=3.52 mm)
- still pending cutting test to 100 mm Y-height

Old FCAL Tungsten Support



- for TB setup nominal combs shifted towards ECAL-P center
- new combs for OLD FCAL tungsten added in the middle (14 slots)
- more mounting holes on the new ECAL-P bottom plate (₱.0 → 1.2mm spacing)


Old FCAL Tungsten Support (zoom)



- for TB setup nominal combs shifted towards ECAL-P center
- new combs for OLD FCAL tungsten added in the middle (14 slots for backup)
- longer ribs no extra side support needed


Old FCAL Tungsten Support

G. Grzelak (University of Warsaw)



- for TB setup nominal combs shifted towards ECAL-P center
- ullet OLD tungsten shifted by 11 mm left o Si sensor coverage
- NOTE: on the drawing all old tungsten plates are cut to 100 mm height

Old FCAL Tungsten (14 plates) plus nominal plates (6)




for TB setup nominal combs shifted towards ECAL-P center

G. Grzelak (University of Warsaw)


ullet OLD tungsten shifted by 11 mm left o Si sensor coverage

NOTE: on the drawing all old tungsten plates are cut to 100 mm height




- BACKUP Solution: instrumented planes with old tungsten
- requires cutting more old tungsten plates to 100 mm Y-height

Position of Silicon Sensors



- instrumented part of the TB setup
- ullet OLD tungsten shifted by 11 mm left (by 2 cells) o symmetric Si sensor coverage
- OLD tungsten: 14cm width : two Si sensors : 18cm width=(only=20-columns for=readout)

Channel map for TB 2025: expected nominal beam aperture



- rotation axis close to geometrical center (decoupling of rotation and translation)
- rotation axis (||) at first instrumented tungsten plate
- readout area: 20 columns × 5.5mm = 110 mm, (15mm OLD tungsten margin on both sides) → a ○


G. Grzelak (University of Warsaw) ECAL-P LUXE 10-APR-2025 1





- Additional plate to host ECAL-P on rotary table
- \bullet dimensions bigger by \sim 10 cm then ECAL-P footprint
- rotation axis at first instrumented plate (close to gravity center of the system)


Support plate + racks scaffolding



- Support plate with ECAL-P and racks scaffolding
- height of the scaffolding "tower" to be adjusted (according to cables length)
- TBD: Faraday Cage and cables routing (light tight throughput) + beam window.

G. Grzelak (University of Warsaw) ECAL-P LUXE 10-APR-2025 13 / 18

Cables: USB-C readout, LV, HV, ... (old CSIS positions!)




TBD: Faraday Cage and cables routing (light tight throughput) + beam window

ECAL-P LUXE

still TBD...

Support plate + racks scaffolding (bottom view)



- Support plate with ECAL-P and racks scaffolding
- additional insulating POM (plastic) plane at the bottom is also needed
- size still under discussion : depends on TB stage table size < ♣ > ♣ < < < <

TB stage at TB22 area (Feb'25)



- ullet table size 120 imes 80 cm, longer side along beam direction
- our "scaffolding" may clash with stage column during rotation...

TB stage: preferred version



- table size 120 × 80 cm (??)
- different table orientation, more space for our setup
- under negotiation with TB crew to have it in TB22 area in June 25 (?)

TB stage: preferred version



- QUESTION: how the X-direction scan was performed?
- can we rely on TB stage translation or should we shift our setup on table ?