## High time resolution & robust longitudinal diagnostics for CLIC

S.P. Jamison, STFC Daresbury Laboratory

W.A. Gillespie, D.A. Walsh University of Dundee









S.P. Jamison / Longitudinal diagnostics WS, DESY Nov 11, 2013

- Electro-optic Transposition (with  $\chi^{(2)}$  interaction)
  - Concept & motivation
  - Experimental conformation of sub systems
- An alignment & tolerance effect in EO systems
- EO-Transposition demonstration & benchmarking @ FLASH
- TFISH THz Field Induced Second Harmonic generation
  - $\chi^{(3)}$  analogue to EO detection
  - Surface effect on Au no phase matching but weak!



## Electro-optic longitudinal diagnostics for CLIC

### Existing systems



- Optical capability limited to >50fs
- Nonlinear/material limited to >60 fs
- "Experimental" in implementation

### EO Optical Transposition

'industrial' nanosecond laser probe

"FROG" for self-referenced <10fs rms

All optical amplification of signal

optical characterisation

Material developments

femtosecond diagnostics without a fs laser

- Target of ~20 fs rms resolution
- Non-destructive
- Instrumentation capability





## Electro-optic longitudinal diagnostics for CLIC

#### Coulomb spectrum shifted to optical region

$$\tilde{E}_{\rm out}^{\rm opt}(\omega) = \tilde{E}_{\rm in}^{\rm opt}(\omega) + i\omega a \tilde{E}_{\rm in}^{\rm opt}(\omega) * \left[\tilde{E}^{\rm Coul}(\omega)\tilde{R}(\omega)\right]$$

### Coulomb pulse replicated in optical pulse

$$E_{\text{out}}^{\text{opt}}(t) = E_{\text{in}}^{\text{opt}}(t) + a \left[ E^{\text{Coul}}(t) * R(t) \right] \frac{\mathrm{d}}{\mathrm{d}t} E_{\text{in}}^{\text{opt}}(t)$$

Monochromatic probe  $\tilde{E}^{\mathrm{opt}}(\omega) \to \tilde{E}_0 \delta(\omega = \omega_0)$ 

$$\tilde{E}^{\text{opt}}(\omega) = \tilde{E}_0 \tilde{E}^{\text{Coul}}(\omega - \omega_0) \tilde{R}(\omega - \omega_0) + \tilde{E}_0 \delta(\omega_0)$$

$$E^{\text{opt}}(t) \sim \left\{ E^{\text{Coul}}(t) * R(t) \right\} \cos(\omega_0 t) + \sin(\omega_0 t)$$

#### Measure spectrum of optical sidebands

 $\left|\tilde{E}^{\mathrm{THz}}(\omega)\right|^2$ 

- Positive and negative sidebands must be same  $\tilde{E}^{\mathrm{THz}}(-\omega) = \tilde{E}^{\mathrm{THz}*}(\omega)$ (some small changes in material response)
- Provides option for intermediate diagnostic (feedback ?)
- test of optics for temporal system & response function

Measure optical pulse envelope in time  $\left\{ E^{\text{Coul}}(t) * R(t) \right\}^2$ 

- FROG diagnostic allows field envelope retrieval  $E^{\text{Coul}}(t) * R(t)$
- Phase of Coulomb converted to envelope of optical bypass FROG carrier-envelope ambiguity



### Final system : single laser delivering ~10ns probe and ~10ns OPA pump



### For testing and specification development



- Ti:S excited THz source
- spectrally filtered TiS for probe,
- synchronized 50ps YAG for OPA pump

## Laser-lab test-bed





S.P. Jamison / Longitudinal diagnostics WS, DESY Nov 11, 2013

## Signal levels, measurability & scaling



Property	Factor of improvement Darsbury -> CLIC
Power <sub>probe</sub>	x36
l	÷100 <sup>2</sup>
r	÷2 <sup>2</sup>
Efield	x186 <sup>2</sup>
Overall	x31

Pulse energy of ~15nJ is produced 1µJ required for single-shot FROG

pulse needs amplifying ~100x

An achievable goal!

### Alignment sensitivity & asymmetry in sidebands





Angular separation of sidebands



- explanation of signal stability / reliability in many previous EO implementations ??
- Mechanism now understood; mitigations through optical design



## Sideband asymmetry in the time domain

Symmetry intrinsic to temporal envelope interpretation *applies to all EO systems, not just narrowband probes* 



Symmetric sidebands spectrum from THz modulation Modified spectrum ; phase and amplitude for causality

Time profile (squared) of modulation Symmetric sidebands spectrum from THz modulation Modified spectrum ; phase and amplitude for causality



## Sideband asymmetry – phase matching with misalignments

For non-collinear THz and probe, positive and negative sidebands diverge

$$\eta_{Intensity} \propto \exp\left(-\frac{1}{2}w^2\left(\Delta k_x^2 + \Delta k_y^2\right)\right) \operatorname{sinc}^2\left(\frac{\Delta k_z L}{2}\right)$$
  
Dominant term for angle

#### Angles involved are small but

- Measurable (<0.1deg for ZnTe)
- potentially distorting on EO
- Dependent on transverse probe dimensions

#### Quantitative characterisation

- imaging spectrometer for sideband angle
- THz angle-only variation







S.P. Jamison / Longitudinal diagnostics WS, DESY Nov 11, 2013

#### **Experimental Data** 4mm ZnTe 0 degrees thz internal angle, 0.5 mm beam width 0 degrees thz internal angle, 0.5 mm beam width 2.5 1.5 0° -1.5 2.5 3--0.1 -0.08 -0.08 -0.04 -0.02 0 0.02 0.04 0.08 0.01 -0.1 -0.08 -0.08 -0.04 -0.02 6 0.02 0.04 0.08 0.08 0.1 0.25 -0.2 -0.15 -0.1 -0.05 ó 1.8871411 degrees thz internal angle, 0.5 mm beam width 3° -1.5 2.5

2.5 1.5 -0.5

-0.08 -0.08 -0.04 -0.02

3-0.1

6°



0 0.02 0.04 0.05 0.08 0.1



**PM Efficiency** 

3.7742821 degrees thz internal angle, 0.5 mm beam width



-0.08 -0.08 -0.04 -0.02 0 0.02 0.04 0.08 0.08 0.1

#### **PM Efficiency** 0.1mm ZnTe change in scales!



1.8871411 degrees thz internal angle, 0.5 mm beam width



3.7742821 degrees thz internal angle, 0.5 mm beam width



THz frequency



THz Angle of incidence

S.P. Jamison / Longitudinal diagnostics WS, DESY Nov 11, 2013



PM Eff. linear fit gradient =  $30.26 \pm 0.04$  THz/° Exp. Data linear fit gradient =  $28.08 \pm 0.10$  THz/°

Discrepancy attributed to potential errors in the positioning of the 20cm angle-position mapping lens and in the position of camera.



## Amplification of post-interaction probe

- Amplification required to allow for FROG characterisation of carved pulse
- NCOPA design demonstrated. 10<sup>3</sup> gain, and sub<20fs rms bandwidth confirmed. (using existing 60 ps YAG pump laser with <ps synchronisation to e-bunch mimic</p>





## System integration

Commercial 'one-package' laser on order - late January 2014 delivery (Litron)

- 10mJ, 10ns , seeded single mode 1062nm + 532nm
- integrated OPA for 750-850nm probe
- Probe bandwidth < 30 GHz</li>
- Stretcher/compressor design in-hand

Modular layout for experimental location flexibility

double layer for compactness? Needed?



### Upconversion spectrogram – a direct X-FROG measurement?



# FROG - spectrally resolved autocorrelation

- *self-referenced* retrieval of temporal profile
- retrieves *envelope* phase within the pulse
- ambiguity remain on carrier-envelope phase. Significant for single-cycle optical pulses
- THz-optical X-FROG: *THz carrier phase ambiguous* 
  - upconversion-FROG: THz carrier phase converted to optical envelope – ambiguity removed



## Material bandwidth limits: Spectral Compositing of Multiple Crystals

*Physics* : *Frequency mixing between Coulomb field (or CSR, CTR,FEL ...) pulse and probe laser* 

Phase propagation mismatch limits bandwidth of materials





#### Material frequency response

Missing response due to phonon resonance High frequency response limited to thin crystals & small signal

- Laser simplicity allows for multiple probes & crystals
- Amplification allows for thin, <50um (?) active thickness crystals
- FROG retrieval of phase information allows for "splicing" of response function



#### Response function for standard crystals

20





• Measure  $\left\{ E^{\text{Coul}}(t) * R(t) \right\}^2$ 

Requires amplification

- FROG retrieval of  $E^{\text{Coul}}(t) * R(t)$
- Evaluate and splice  $E^{\tilde{Coul}}(\omega)$

Currently do not have very short test pulses - will demonstrate FROG on single xtal output

+ GaSe for shorter pulses; BBO in extreme case (? on validity in extreme)



S.P. Jamison / Longitudinal diagnostics WS, DESY Nov 11, 2013

## Spectral Transposition

Ideally, nonlinear  $\chi^{(2)}$  frequency mixing transposes complex amplitude to the optical regime.



Optical measurement techniques such as **FROG** can be used to retrieve full spectral **amplitude** and **phase**.

However, material response induces distortions and limits bandwidth



## Material Limits

#### 100 $\mu m$ ZnTe and GaP





#### 10 $\mu m$ ZnTe and GaP





S.P. Jamison / Longitudinal diagnostics WS, DESY Nov 11, 2013

## Higher bandwidth materials - TFISH THz Field Induced Second Harmonic Generation

Gold mirror

- THz characterisation on metal surface demonstrated e.g. Nataha (1999)
  - limited examples, efficiency needs characterisation
  - will need (extra) amplification
- separate ASTeC project in gold surface DFG for UV characterisation

- being adapted for TFISH tests with laser generated THz

- Enhancement of interaction: known enhancement of SHG with nanoparticle metals.
  - Dundee producing samples, Daresbury characterisation Marie Curie student (Mateusz Tyrk);
  - Silver-in-glass MDNs generated by field-assisted ion exchange (30nm dia nanoclusters)
  - Nanoparticles re-shaped to ellipsoidal shape using polarised ps laser fluence







From Nataha et al.

Spectronice Metal SHG set-up - Being adapted for CLIC TFISH

## Accelerator tests, and benchmarking with short pulses



## Laser lab tests limited by THz source

- Bandwidth <6THz
- Calibrated spectrum & temporal profile missing characterised by EO; circular "benchmarking

## CTR experiments at ALICE (unsuccessful)

- no CTR transport system,
- without amplification
- v. low CTR pulse energy (25MeV, 60pC)

CTR pulse carving experiments (ALICE, Daresbury lab.)

## **Proposal:** Benchmarking against FLASH transported CTR

- Cross check against calibrated CTR spectrometer
  - Individual crystals for selected bandwidth regions
  - Direct spectrum comparison
- Evaluate temporal retrieval with realistic field
  - Compare with expected bunch profile & CTR propagation?
  - *Icing*: compare with 'simultaneous' TDC measurements



**Proposal:** Benchmarking against FLASH transported CTR

### Optical arrangement

- Compact layout ~ 1.3x1.0m (stacking)
- Open layout 2.4 x 1.0 m
- Optical window access for probe to CTR
- Focal spot for CTR (within CTR spectrometer?)
- Optical flipping between spectrometer / EO measurements

### Installation and build-up

- Daresbury/Dundee provide system on breadboard
- David Walsh available for extended periods setup and experiment
- Additional Daresbury people available for 1-2 week experiment runs
- CTR spectrometer and TDC for experiments
- Laser expected late January 2014; system build by end March (linked to laser lab relocation schedule)
- Available for DESY/FLASH experiments ~June 2014



### Optical arrangement – "open layout"



#### + air cooled chillier



### In conclusion

- Electro-optic transposition scheme being developed for CLIC
  - fs resolution without fs laser
  - Concept allows multiple crystal detectors
  - Component systems tested
- Integrated 'one-laser' system ready ~march 2014
- Separate project on TFISH surface effects
- On-going materials development at Dundee, non-linear testing at Daresbury started

# Thank you.



Experimental spectrograms showing angular dispersion of sidebands at several THz input angles of incidence



-0.5

-1.5

-2 -

-2.5 -

0.1 0.08 0.06 0.04 0.02 0 0.02 0.04 0.06 0.08 0.1

-0.5

.15

-2

-2.5 -

0.1 0.08 0.06 0.04 0.02 0 0.02 0.04 0.06 0.08 0.1

-0.5

.15

-2

-2.5

0.1 0.08 0.06 0.04 0.02 0 0.02 0.04 0.06 0.08 0.1

.0 5

.15

-2.5

0.1 0.08 0.06 0.04 0.02 0 0.02 0.04 0.06 0.08 0.1

## Frequency resolved optical gating (FROG)

#### Ultrafast laser technique, self referenced characterisation of pulse at <10fs res.

Upconverted optical pulse is too weak for FROG few nJ available, uJ required. Solution: Non-collinear Chirped Pulse Amplification (NCPA)



#### Commercial packaged laser specification

